- Ⓡ FR APPLICATEURS DE CLIPS HEMOSTATIQUES
- Ⓡ EN HEMOSTATIC CLIP APPLIERS
- Ⓡ CZ APLIKÁTOR HEMOSTATICKÝCH SVOREK
- Ⓡ DA HÆMOSTATISKE KLEMMER
- Ⓡ NL APPLICATORS VOOR HEMOSTATISCHE CLIPS
- Ⓡ DE HÄMOSTATISCHE CLIPAPPLIKATOREN
- Ⓡ EL ΛΑΒΙΔΑ ΑΙΜΟΣΤΑΤΙΚΩΝ ΚΛΙΠ
- Ⓡ HU VÉRZÉSCSILLAPÍTÓKLIP-APPLIKÁTOROK
- Ⓡ IT APPLICATORI DI CLIP EMOSTATICHE
- Ⓡ LT HEMOSTATINIŲ SPAUSTUKŲ APLIKATORIAI
- Ⓡ NO HEMOSTATISKE KLIPSAPPLIKATORER
- Ⓡ PL APLIKATORY KLIPSÓW HEMOSTATYCZNYCH
- Ⓡ PT PINÇAS HEMOSTÁTICAS
- Ⓡ RO APLICATOARE DE CLEME HEMOSTATICE
- Ⓡ SK APLIKÁTORY HEMOSTATICKÝCH SVORIEK
- Ⓡ ES APLICADORES DE CLIPS HEMOSTÁTICOS
- Ⓡ SE HEMOSTATISKA TÄNGER
- Ⓡ AR أدوات مسك الملاقط الجراحية المُرَقنة
- Ⓡ ZH 止血夾施夾器
- Ⓡ RU ГЕМОСТАТИЧЕСКИЕ КЛИП-АППЛИКАТОРЫ
- Ⓡ TR HEMOSTATİK KLIPS APLIKATÖRÜ

FRANÇAIS .....	3
ENGLISH .....	6
ČEŠTINA .....	9
DANSK .....	12
NEDERLANDS .....	15
DEUTSCH .....	18
ELAHNIKA .....	21
MAGYAR .....	24
ITALIANO .....	27
LIETUVIŲ K .....	30
NORSK .....	33
POLSKI .....	36
PORTUGUÊS .....	39
ROMÂNĂ .....	42
SLOVENČINA .....	45
ESPAÑOL .....	48
SVENSKA .....	51
العربية .....	54
繁體中文 .....	57
РУССКИЙ .....	60
TÜRKÇE .....	63
SYMBOLS .....	66

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# HEMOSTATIC CLIP APPLIERS

## DEVICE IDENTIFICATION

**Hemostatic clip applicators** are Class 1r stainless steel devices. They are reusable non-sterile surgical instruments used in open surgery. A Unique Identifier (UDI) in the form of a Datamatrix code and a series of numbers are affixed to the instrument in the central part.

## INTENDED USE

**Hemostatic clip applicators should only be used with VITALITEC SLS-CLIP®, VITALITEC CLIP 9 and CLIP GEM™ in open surgery.**

The 7 available sheath colours allow the identification of hemostatic clip applicators adapted to the size of the **VITALITEC SLS-CLIP®, VITALITEC® CLIP 9 and CLIP GEM™**.

## INDICATIONS

**Hemostatic clip applicators** are designed to ensure optimal grip, hold and closure of **VITALITEC SLS-CLIP®, VITALITEC® CLIP 9 and CLIP GEM™**.

## TARGET POPULATION

**Hemostatic clip applicators** must only be used by surgeons. Users must be trained in surgical techniques and procedures involving the combined use of Hemostatic clip applicators and **VITALITEC SLS-CLIP®, VITALITEC® CLIP 9 and CLIP GEM™**.

**Hemostatic Clip Applicators** are intended to be used in all patients regardless of age or sex.

## COMPOSITION

**Hemostatic clip applicators** are made of stainless steel; the rings are coated with a coloured polyamide sheath.

## MODE OF ACTION

The mechanical compression exerted by the Hemostatic clip applicators allows the **VITALITEC SLS-CLIP®, VITALITEC® CLIP 9 and CLIP GEM™** hemostatic clips to press biological tissue.

## CLINICAL PERFORMANCE AND CLINICAL BENEFITS

The clinical benefits and clinical performances of the **Hemostatic clip applicators** depend on the **VITALITEC SLS-CLIP®, VITALITEC® CLIP 9 and VITALITEC® CLIP GEM™ hemostatic clips**.

## INSTRUCTIONS BEFORE EACH USING

After cleaning and before sterilization, it is recommended to treat the applicators and especially the connection areas as well as all moving parts with a biocompatible lubricant (silicone-free, oil-free and water-based) suitable for surgical instruments before sterilization.

**The service life of hemostatic clip applicators** is not defined in time or limited by a number of reprocessing cycles. Like any surgical instrument, **hemostatic clip applicators** are subject to wear and tear from normal use. It is therefore imperative that each instrument be thoroughly inspected to ensure that it is functioning properly. Each hemostatic clip applicator should be visually inspected before each use, checking the lack for:

- Dirt on all surfaces, joints and moving parts.
- Signs of abnormal wear, cracking or corrosion.
- Rust due to the appearance of microscopic cracks.

And by checking also the state:

- Moving parts, pins, springs, grooves, connection areas.
- No flaking of the colored ring sheath coating.

Each **hemostatic clip applicator** should be functionally inspected to ensure proper operation before each use, checking for:


- Slack in the moving parts.
- Broken or missing parts.

Particular attention must be paid to the jaws of the instrument. If an instrument shows any abnormality, **BE CAREFUL, DO NOT USE IT**

- Make sure that the two jaws, when closing, first touch each other distally (A) and are not twisted (B)

- Ensure that the jaws are correctly aligned to avoid shearing or poor closure of the **VITALITEC SLS-CLIP®, VITALITEC® CLIP 9 and GEM™ CLIP (C)**.

- Make sure that the grooves of the jaws are not deformed (D) and/or do not have any impact (E).

## INSTRUCTIONS CLEANING AND STERILIZATION

Each **hemostatic clip applicator** should be cleaned and sterilized before use. The plastic jaw protection cap and primary packaging must be removed prior to the first cleaning.

### INFORMATIONS ON VALIDATED CLEANING PARAMETERS:

The manual and automatic cleaning parameters have been validated with Neodisher Medizym enzymatic detergent.

The final rinse is performed with deionized water with a low bacterial load.

## GENERAL REMARKS

This product is delivered and packaged non-sterile, it must be cleaned and sterilized before use.

- The brushes are involved in the effectiveness of the cleaning. They must be inspected and replaced, if necessary. Do not use metal brushes. Use soft bristles.
- The cleaning equipment must be validated according to ISO15883 and sterilization equipment according to ISO17665 standard and be subject to regular maintenance and controls. These validations are the responsibility of the users.
- For cleaning agents, make sure to follow the manufacturer's instructions regarding concentration, contact time and temperature.
- For water, use running water for intermediate rinses and deionized water for final rinses,

## PERSONNEL PROTECTION:

Personnel must wear gloves when handling the cleaning products, (Refer to the safety data sheets of the products used).

## PREPARATION OF THE INSTRUMENT:


Before cleaning and sterilization, the **hemostatic clip applicator** must be opened manually. To do this:

- Disengage the center spring.
- Spread the legs to the maximum extent.

The spring must be disengaged by hand, not with an instrument.

In order to avoid its deformation that could affect the settings of the jaws.

It is important not to stack the instruments on top of each other instruments on top of each other to avoid not to damage the jaws.

## PRE-CLEANING:

The pre-decontamination operation should be performed within 4 hours of using hemostatic clip applicators:

- 1- Open the applicator.
- Remove major impurities immediately after use using a disposable cloth.
- 2- Immerse the applicator for 10 minutes in an enzymatic detergent or equivalent according to the temperature, concentration and quality of the water recommended by the manufacturer.
- 3- Rinse thoroughly (1 minute minimum) under tap water.
- 4- Place the device in an ultrasonic bath filled with fresh enzymatic detergent or equivalent prepared according to the manufacturer's instructions for water concentration and quality. Sonicate for a minimum of 10 to 15 minutes at a frequency of 35 to 45 kHz, the temperature should be 20°C (not exceeding 40°C).
- 5- Rinse thoroughly (minimum 1 minute) with tap water.

## MANUAL CLEANING:

- 1- Immerse the applicator for at least 15 minutes in an enzymatic detergent or equivalent, prepared according to the manufacturer's instructions for temperature, concentration and water quality. While immersed, brush the entire applicator with a soft bristle brush for 1 minute (minimum) or until all visible soil is removed.
- 2- Rinse thoroughly with tap water and scrub with a clean, soft bristle brush for at least 2 minutes.
- 3- Immerse in a bath of clean running water and let soak for at least 2 minutes.
- 4- Rinse thoroughly (1 minute minimum) under deionized water
- 5- Dry the instruments with a clean, lint-free cloth.
- 6- Check that there is no visible dirt.

## AUTOMATIC CLEANING:

Use an enzymatic detergent or equivalent suitable for the washer-disinfector, following the manufacturer's recommendations for concentration, time and temperature. The recommended automatic cleaning parameters below have been validated with Neodisher MediZym detergent. Rinse water must be deionized.

- 1- Place the applicator in a suitable holder; instruments should not be stacked or in contact with each other.
- 2- Start a Vario TD or equivalent program cycle according to the recommendations below.
- 3- Remove the applicator from the washer-disinfector

Stage	Recirculation time	Temperature	Detergent Type
Pre-wash	3 min	Cold Water	Not applicable
Cleaning	10 min	50°C water*	Enzymatic detergent or equivalent
Rinsing	6 min	Cold Water	Not applicable
Disinfection thermal	5 min	90 °C	Not applicable
Drying	35 min	100°C	Not applicable

\* Refer to manufacturer recommendations for optimal temperature. Do not exceed 50°C.

It is important to lay the **hemostatic clip applicators** flat, with the jaws facing up, to prevent them from getting stuck in the basket of the automatic washer.

## CONTROL:

After cleaning-disinfection operations, check the applicator for corrosion, dirt or damaged parts. Special attention should be paid to the applicator jaws.

**When lubricating the instrument**, use a water-soluble medical lubricant suitable for steam sterilization and surgical instruments.

## PACKAGING:

Place the product in a single package in a sealed pouch suitable for steam sterilization meeting ISO 11607-1 or applicable local standards.

## STERILIZATION:

Use a vacuum sterilizer in accordance with ISO17665.  
Steam sterilization to be implemented according to the following cycles:  
Sterilization cycle: 132°C for at least 4 minutes or 134°C for at least 3 minutes (SAL 10-6). **France**, validated cycle: 134°C for 18 minutes.  
Drying time: 30 minutes minimum.

## PRECAUTIONS FOR USE

### CHECK BEFORE USE:

Visual and functional checks before using the instrument must be carried out to ensure its correct operation.

### CLOSING FORCE:

**Be sure to apply firm pressure on the rings to overcome** the spring force. Insufficient pressure may result in incomplete closure of the clip.

### USE FOR INTENDED PURPOSE:

Improper use of the instruments may result in stress and irreversible damage.

### HANDLE WITH CARE:

Improper handling may result in contamination, corrosion, and cracking.  
Clip appliers should be handled with extreme care to prevent the clip from falling out of the jaws.  
The jaws should be inserted vertically until they stop at the bottom of the cartridge so that the clip is loaded and in the correct position in the grooves.

### AVOID DRIED RESIDUE:

Blood and other secretions can cause discoloration and corrosion.

### AVOID LEAVING INSTRUMENTS WET:

Instruments left in liquid can cause a galvanic current that can lead to corrosion.

### AVOID MECHANICAL STRESS:

Forcing the opening to a maximum amplitude and/or applying too much pressure at closing can cause misalignment of jaws as well as stress corrosion.

### FOLLOW PROCESSES:

Ensure that sterilization department routines are strictly followed.

## CONTRAINDICATIONS

**Hemostatic clip appliers** are not intended for any use other than that indicated.  
In the case of a patient with a prion disease, the device should not be reused under any circumstances.

## UNDESIRABLE EFFECTS

No known undesirable effect related to use of Clips Appliers is reported. Refer to the instructions for uses of the **VITALITEC SLS-CLIP®, VITALITEC® CLIP 9** and **VITALITEC® CLIP GEM™** hemostatic clips for adverse events to use of clips.

Any serious incident related to the use of this device should be reported to the competent authority of the country in which the user and/or patient is located and to the manufacturer or his representative.

## STORAGE/HANDLING/ELIMINATION

When stored, this device should be kept in a dry, dark place and under conditions that maintain the functional properties of the device.

Never attempt to repair an instrument. In the event of instrument failure, please contact PETERS SURGICAL or its local representative to determine the appropriate action to be taken.

Any used device should be disposed of in an appropriate incinerator in accordance with local or national regulations for the disposal of hospital waste.

## INFORMATION TO BE PROVIDED TO THE PATIENT

The applier is an accessory to the Clips.  
Please refer to the instructions for the clips for the transmission of information to the patient by the surgeon concerning the contraindications, side effects and risks associated with the use of this device.

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