ENGLISH		BASIC UDI-DI		BASIC UDI-DI	
CYGNET® Instrument Tray		085798706-90046UA		085798706-90046UA	
Rx Only	Prescription Only	UDI	Unique Device Identifier	UDI	Unique Device Identifier
Non Sterile	Distributor	MD	Medical Device	MD	Medical Device
Manufacturer	REF	Catalog Number	Consult Instructions for Use	Manufacturer	REF
Country of Manufacture	LOT	Batch Code		Country of Manufacture	LOT
EC REP	Authorized representative in the European Community	CH REP	Authorized representative in Switzerland	EC REP	Authorized representative in the European Community

INDICATIONS FOR USE

The Cygnet® surgical instrument tray is to be used for sterilizing and storing the Cygnet® surgical instruments. The Cygnet® surgical instrument tray must be used in conjunction with a sterilization wrap that is cleared by the regulatory authorities for the indicated sterilization cycle.

INTENDED PURPOSE

The Cygnet® surgical instrument tray is intended to enclose and protect the Cygnet® clamp during sterilization and storage.

CONTRAINDICATIONS

WARNINGS

- The Cygnet® surgical instrument tray and mat are provided NON-STERILE and must be cleaned and sterilized according to these instructions prior to use.
- Stacking of trays will adversely affect sterilization and drying effectiveness. DO NOT STACK trays in the sterilization chamber.
- Do not put more than one Cygnet® clamp in a tray.
- Ethylene oxide (EO), gas plasma and dry heat sterilization methods are not recommended for sterilization of Cygnet® surgical instrument tray. Steam (moist heat) is the recommended method.
- Personal Protective Equipment (PPE) should be worn when handling or working with contaminated instruments or trays.
- Saline and cleaning/disinfection agents containing aldehyde, chloride, active chlorine, bromine, iodine or iodide are corrosive and should not be used.
- Use of hard water should be avoided. Softened tap water may be used for most rinsing however purified water should be used for final rinsing to prevent mineral deposits.
- Proper disposal of contaminated devices is required, including disposal of reusable instruments that are no longer acceptable for use.

PRECAUTIONS

INTENDED USERS

Professionals trained in surgical instrument cleaning, disinfection and sterilization

PATIENT TARGET GROUP

Pediatric to Geriatric

COMPLICATIONS/ADVERSE EVENTS

Cross contamination leading to infection

HOW SUPPLIED

The Cygnet® surgical instrument tray is supplied non-sterile, one per package. This can be ordered with or without the Cygnet® clamp.

INSTRUCTIONS FOR USE

These instructions are in accordance with ISO 17664 and AAMI S781. They apply to Cygnet® surgical instrument trays supplied by Vitalizec International and intended for reprocessing in a health care facility setting.

Cleaning/Sterilization

CAUTION: Vitalizec recommends the use of enzymatic detergents for optimum cleaning results. Vitalizec also anbefaler at bruge enzymatiske rengøringsmidler for at opnå optimale rengøringsresultater.

CAUTION: Trays are reusable but may be damaged after multiple uses. Always inspect tray components for cracking, chipping or other signs of damage before use. **CAUTION:** Only use accessory components that have been designed for the Cygnet® surgical instrument tray.


CAUTION: Do not use solvents, abrasive cleaners, metal brushes or abrasive pads when cleaning the Cygnet® surgical instrument tray.

Manual Cleaning

- Disassemble components prior to cleaning.
- Put the tray in cold running utility (tap) water to remove gross soil.
- Prepare a pre-soaking, neutral pH enzymatic cleaner bath such as Enzol at the manufacturer's recommendation of 1 oz per gallon.
- Soak the tray for at least 10 minutes. While soaking, scrub surfaces using a soft bristled brush until all visible soil has been removed. Pay particular attention to cleaning between the pins on the mat.
- Carefully inspect and repeat cleaning sequence if any visual debris remains on any of the exposed surfaces.
- Dry the surfaces of the tray with a clean, lint free cloth.

NOTE: The Cygnet® surgical instrument trays have been validated for use in a Pre-Vacuum steam sterilization environment. The validation testing was performed by an independent and fully accredited testing laboratory. For additional information or to request a copy of the report please contact Vitalizec International.

Configure the tray and mat properly to accommodate the Cygnet® clamp to be sterilized. The mat should be used to help stabilize and protect the clamp. See Figure 1 below.

- Wrap the Cygnet® surgical instrument tray with an approved non-woven textile wrap. The trays should always be placed flat on the sterilizer rack prior to sterilization.
- Stacking of the trays inside the autoclave chamber is not recommended.
- Sterilize the load. Validated exposure times and temperatures to achieve a 10⁶ sterility assurance level (SAL) are listed below.

Pre-vacuum Cycle:

Preconditioning Pulses:	4
Temperature:	132° C / 270° F
Full Cycle Time:	4 Minutes
Dry Time:	30 Minutes

- A 30-minute cooling time is recommended after drying but longer times may be necessary because of load configuration and ambient temperature and humidity.
- If wet spots on the packaging are observed, the contents must be treated as non-sterile and be thoroughly reprocessed.

Preconditioning and drying must be validated by the institution performing the cleaning and sterilization processes.

Inspection
Always inspect tray components for cleanliness and confirm that there is no visual contamination prior to use.

STORAGE AND HANDLING
The Cygnet® surgical instrument tray should be stored in an area that is dry and clean, away from heat and chemical fumes.

CLINICAL BENEFITS
There are no specific clinical benefits for the instrument tray. The tray has no clinical application. It is a device used only to store and sterilize the Cygnet Clamp.

NOTICE

If any serious injury has occurred to the patient or user of this device, notify Vitalizec and the competent authority of the Member State in which the user and/or patient is established.

PRODUCT WARRANTY
All Vitalizec Cygnet® surgical instrument trays are provided "AS-IS". Customer is responsible for any repair or replacement.

PRODUCT INFORMATION DISCLOSURE
Vitalizec excludes all warranties, whether expressed or implied by operation of law or otherwise, including, but not limited to: Any implied warranties of MERCHANTABILITY or FITNESS, since storage and handling of this device by the user, as well as factors relating to the patient, the surgical technique, the surgical therapy, and other factors beyond Vitalizec's control directly affect this device and the results obtained from its use. Vitalizec shall not be liable for INCIDENTAL or CONSEQUENTIAL LOSS, DAMAGE, or EXPENSE directly or indirectly arising from the use of this device.

Vitalizec neither assumes, nor authorizes any other person to assume for it, any other form of LIABILITY or RESPONSIBILITY in connection with this device.

	EC REP
VITALITEC INTERNATIONAL, INC. d/b/a Peters Surgical 10 Cordage Park Circle Suite 100/200 Plymouth, MA 02360 USA Phone: 508-747-6033 Customer Service: 800-805-6590 www.peters-surgical.com	EMERGEO EUROPE Westvortsevestdijk 60, 6827 AT Arnhem The Netherlands

	CH REP
DISTRIBUTED BY: PETERS SURGICAL immediate AUREUM 1 cours de l'le Seguin 92100 Boulogne - Billancourt TEL: 33 (0) 1 48 10 6262 commercial@peters-surgical.com	MEDENOVIO SWITZERLAND Gothardstrasse 28 6302 Zug Switzerland

CLINICAL BENEFITS
There are no specific clinical benefits for the instrument tray. The tray has no clinical application. It is a device used only to store and sterilize the Cygnet Clamp.

NOTICE
If any serious injury has occurred to the patient or user of this device, notify Vitalizec and the competent authority of the Member State in which the user and/or patient is established.

PRODUCT WARRANTY
All Vitalizec Cygnet® surgical instrument trays are provided "AS-IS". Customer is responsible for any repair or replacement.

PRODUCT INFORMATION DISCLOSURE
Vitalizec excludes all warranties, whether expressed or implied by operation of law or otherwise, including, but not limited to: Any implied warranties of MERCHANTABILITY or FITNESS, since storage and handling of this device by the user, as well as factors relating to the patient, the surgical technique, the surgical therapy, and other factors beyond Vitalizec's control directly affect this device and the results obtained from its use. Vitalizec shall not be liable for INCIDENTAL or CONSEQUENTIAL LOSS, DAMAGE, or EXPENSE directly or indirectly arising from the use of this device.

Vitalizec neither assumes, nor authorizes any other person to assume for it, any other form of LIABILITY or RESPONSIBILITY in connection with this device.

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