
Peters SURGICAL
 Immeuble AURELIUM
 1 Cours de l'Île Seguin
 92100 BOULOGNE-BILLANCOURT
 FRANCE
 Tél. +33 (0)1 48 10 62 62
 Fax +33 (0)1 48 91 22 99
 E-mail : peters@peters-surgical.com
www.peters-surgical.com

MD REF MAUT60

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INSTRUCTIONS FOR USE – PETERS SURGICAL SINGLE-USE STERILE SURGICAL UTERINE MANIPULATOR

CAUTION

Prior to using the PETERS SURGICAL single-use sterile surgical uterine manipulator, please read these Instructions carefully. Failure to comply with the following recommendations may entail serious consequences. Laparoscopic operations should only be performed by qualified doctors having been specifically trained for the purpose and having gained sufficient experience of use. All surgeons performing laparoscopic operations should be aware of the possible risks and complications involved when using this surgical technique (notably risks for the patient, as well as in manipulating instruments connected to an external power source, such as electrosurgery, ultrasound, laser, gas and liquid pressure). Adequate details of the operation should be given to all patients, such as provided in current legislation with regard to the right to information.

DESCRIPTION OF THE MEDICAL DEVICE

List of all elements contained in the sterile packaging:

- 1 dilatation bougie;
- 1 uterine manipulator equipped with:
 - a cutting-guide valve, and
 - a balloon;
 - 1 balloon handle.

The single dilatation bougie unit enables to calibrate the cervix and to assess the length of the latter, as well as that of the uterine cavity. Due to the shape of the bougie, dilatation is conical and is conducted in stages. The surgical uterine manipulator is comprised of a metal tube with a sharp screw tip on one end, and a handle for manipulation on the other. The cutting-guide valve may freely rotate and translate, thereby enabling to ensure precise cutting (removal of the uterus). The balloon enables to guarantee pneumoperitoneum as soon as removal of the uterus begins. At the end of surgery, the balloon handle (linked to the balloon) enables to maintain pneumoperitoneum.

► On no account may PETERS SURGICAL be held liable in the case of failure by the operator to observe proper surgical practices.

SURGICAL UTERINE MANIPULATOR INDICATIONS

The PETERS SURGICAL surgical uterine manipulator is a single-use medical instrument, exclusively designed for gynaecological surgery. It should only be used by qualified practitioners having gained due experience.

The surgical uterine manipulator is designed to expose the uterus (for moving in and around the abdominal cavity). It further enables to facilitate presentation of the rear of the vagina during a full hysterectomy, while preserving pneumoperitoneum.

CONTRAINDICATIONS

Use of the PETERS SURGICAL surgical uterine manipulator is strictly forbidden if there exists a contraindication in relation to laparoscopic surgery.

Do not use the device for pregnant patients, whether assumed or confirmed, nor in the case of a pathology not compatible with a uterine cannulation.

► On no account may PETERS SURGICAL be held liable for assessment of the contraindication in relation to laparoscopic surgery.

WARNINGS

► These instructions solely refer to information concerning the PETERS SURGICAL surgical uterine manipulator and to its method of use. On no account should they be considered as an instructions' manual for using surgical uterine manipulators in general.

► PETERS SURGICAL sterile surgical uterine manipulators are designed for single use only. On no account may they be used a second time over, whatever the reason or the means.

► PETERS SURGICAL surgical uterine manipulators are delivered in sterile packaging. They should never be resterilised, regardless of whether they have been used or not.

► Any reuse or resterilisation of this instrument shall affect its original sterility and may lead to an infection or to cross-contamination in patients. Reuse or resterilisation may also compromise the integrity of the instrument and the use thereof, possibly causing injury, sickness or even death of the patient.

► Warnings specific to a particular manipulation are stated herein, throughout the present Instructions (in bold type and preceded by ►).

► Responsibility with respect to the storage, the opening, the appropriate use, the handling, the insertion and the destruction of the surgical uterine manipulators lies solely with the operator. On no account may PETERS SURGICAL be held liable in the case of failure by the operator to observe proper surgical practices for the storage, the opening, the appropriate use, the manipulation, the insertion and the destruction of the surgical uterine manipulators. Likewise, the manufacturer may not be held liable in the case of reuse and/or resterilisation.

► Never use the sterile packaging with its contents in the following cases (whereupon it should be destroyed):

- Packaging severely crushed or clearly damaged;
- Visible damp marks on the cardboard packaging;
- Perforation of the outer covering.

STORAGE CONDITIONS

- Storage temperature: avoid any prolonged exposure to low and/or high temperatures.
- Store in an environment adapted to sterile medical devices.
- For transportation, all the above points should be applied.

Information for the United States (USA):

- US Federal Law only allows the sale of this instrument by a doctor or upon presentation of a doctor's prescription.
- Sterilised with ethylene oxide.
- Single use only - Not to be reused.

In the event of any unfortunate occurrence (or vigilance on equipment), please contact PETERS SURGICAL at the following address: peters@peters-surgical.com

INSTRUCTIONS FOR USE

Prior to using the PETERS SURGICAL surgical uterine manipulator, the user should ensure that the operating zone, i.e. the abdominal wall, the cervix and the vaginal cavity, is regularly prepared according to the rules for proper surgical practices. Such operations involve inserting the bougie and the manipulator into both the cervix and the uterus via the body's natural channels, by adopting regular gynaecological techniques in order to protect the surrounding tissue from any possible injury.

1. Open the bulkpack casing and withdraw one of the medical instruments in its sterile packaging.

2. Make sure that the packaging guaranteeing sterile conditions for the medical instrument is neither open nor damaged.

3. Check the useby date printed on the packaging.

4. Open the sterile packaging according to the proper practices for the relevant operation zones.

► PETERS SURGICAL surgical uterine manipulators and their packaging should be destroyed as of when sterility may no longer be guaranteed, notably in the following cases: expiry of the useby date on the packaging; packaging severely crushed or clearly damaged; damp marks visible on the packaging (both the cardboard box and the outer covering); perforation of the outer covering.

5. Check:

- the good state of the surgical uterine manipulator;
- the presence of all elements (dilatation bougie, manipulator, cuttingguide valve, balloon and balloon handle).

► Damaged surgical uterine manipulators should not be used.

6. Before starting surgery, make sure that all the laparoscopic instruments are adapted to PETERS SURGICAL surgical uterine manipulators. Check the state and the working order of all instruments prior to use.

► On no account may PETERS SURGICAL be held liable in the case of failure by the operator to observe proper surgical practices when verifying whether the laparoscopic instruments are compatible with PETERS SURGICAL surgical uterine manipulators.

7. Before inserting the bougie and the surgical uterine manipulator, verify the patient's internal body conformation in order to prevent any risk of injury, whether of a vascular nature or to the vital organs.

► Caution: excessive pushing or bending when setting in place or manipulating may damage the organs and tissues (risk of uterine perforation).

8. Position a speculum in order to open up the region of the cervix by observing proper surgical practices.

9. Grasp the cervix using Pozzi forceps by observing proper surgical practices.

10. Insert the narrow end of the PETERS SURGICAL singleuse bougie and start dilating the cervix in successive stages (conical dilatation) by observing proper surgical practices. Repeat the operation using the wide end of the PETERS SURGICAL bougie. This shall enable the correct calibration of the cervix for the subse quein insertion of the manipulator. The graduated PETERS SURGICAL singleuse bougie shall provide an indication of the depth of the uterine cavity.

11. Grasp the handle of the PETERS SURGICAL surgical uterine manipulator.

12. Screw (clockwise) the manipulator's threaded tip into the cervical canal, held in place by Pozzi forceps, thus enabling the cone to lean against the cervix. Then remove the forceps.

13. Insert the cuttingguide valve and the non-inflated balloon into the vagina.

14. The PETERS SURGICAL surgical uterine manipulator is now in position to enable the various stages of hysterectomy dissection under coelioscopy, pursuant to proper surgical practices.

► Electrosurgery apparatuses (or any other energy source, e.g. ultrasound, laser) should be used at minimal power, though under constant direct visual control (until reaching the intraabdominal side of the uterus). Avoid any direct contact from such energy sources with the manipulator, the heat released possibly damaging the valve or the tip.

15. Before beginning to open the vagina, inflate the balloon using a syringe with approximately 20 ml of air, then adapt the air volume to the patient's morphology (see sketch no. 6) in order to guarantee pneumoperitoneum as from the start of incision.

► In the event of expulsion of the balloon, never attempt to replace it before having first deflated it; then reinsert the balloon and inflate it again with the necessary volume of air.

16. Translate and revolve the cuttingguide valve using the graduated handle (mark 7A on the serrated wheel should be in line with mark 7B on the valve: see sketch no. 7) in order to show the vaginal Douglas culde sacs and to enable vaginal incision across the whole circumference (360°).

17. Once the circular incision terminated :

a) if the size of the uterus enables extraction via the vaginal route, remove all the instruments from the trocars, stop insufflations, deflate the balloon and extract the operated part. Unscrew (anticlockwise) the threaded tip from the uterus;

b) if the size of the uterus prevents extraction via the vaginal route, cut it into pieces using regular surgical practices. Deflate the balloon and withdraw the manipulator.

18. Separate the cutting guide valve into two parts using a translation gesture and release the balloon (see sketch nos 1 to 3).

19. Fix the balloon handle onto the balloon (see sketch nos 4 and 5).

20. Insert this balloon assembly into the vagina.

21. Inflate the balloon using a syringe (see sketch no. 6) with approximately 20 ml of air, then adapt the air volume to the patient's morphology in order to enable reconstitution of the pneumoperitoneum and continue the laparoscopy.

► In the event of expulsion of the balloon, never attempt to replace it before having first deflated it; then reinsert the balloon and inflate it again with the necessary volume of air.

22. Suture the rear of the vagina using regular surgical techniques.

23. Once the operation terminated, deflate the balloon and extract the balloon assembly.

24. A gynaecological examination should be systematically performed via vaginal palpation at the end of the operation.

25. Make sure that the surgical uterine manipulator is fully intact, with notably the presence of both parts of the cuttingguide valve, as well as the balloon and the balloon handle.

After surgery, the surgical uterine manipulator and its accessories should be disposed of in the same manner as all operatingtheatre waste.

► The surgical uterine manipulator, its accessories and its packaging should be destroyed after use.

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