**SILISOFT® SILICONE WOUND DRAINS (ROUND OR FLAT)****DESCRIPTION**

Silicone fluted drain round or flat, with or without stainless steel trocar. Connector for drainage bottle added with each drain. X.R.O. line

Round fluted drain: total length around 110 cm. Fluted over around 20 - 30 cm. FG 10 – 15 – 19 – 24.

Flat ¾ fluted drain: total length around 110 cm. Fluted over around 20 cm. Width: 7mm – 10mm

Flat fully fluted drain: total length around 102 cm. Fluted over around 20 cm. Width: 7 mm – 10 mm

**INDICATIONS**

Indicated for use in surgical situations where small amounts of drainage are expected. Surgical drainage under depression.

**ACTIONS**Use in peroperative:

Remove the cap using a rotating motion. The fluted part of the drain is inserted into the wound or the operative cavity. The drain trocar perforates the skin several centimetres from the inside of the wound towards the outside. The mark at 5cm after the end of opened channels facilitates the positioning of the drain in the surgical field. The drain is fixed to the skin by means of a thread without tightening excessively, caution should be taken not to perforate the drain with the needle or any other sharp instrument what would result in damaging the tube. On the assumption that the end would be shortened in per operative, it would be advisable to note it in the patient file in order to take it into account during the ablation of the drain.

The trocar is eliminated after use.

Mobilisation and removing as directed by the surgeon.

SILISOFT® drains must not be used for more than 10 days in continuity

In case of connection with a Redon bottle, proceed as follows:

Connector of Drainage Bottle	Cut at level FG 8 or 10	Cut at level FG 12	
Drain and connector	FG10 connector 2.8 x 2.8 mm	FG15-19 connector 4 x 4 mm	Width 7-10 mm connector 4 x 4 mm.

**CONTRAINDICATIONS**

None contraindication known.

**ADVERSE EVENTS**

None.

**CAUTIONS**

Do not suture through these drains, or handle with any instruments! During placement the drains, as this may lead to breakage. Drains should be placed and removed carefully, by applying a slow steady pressure. Excessive force may cause the drains to break. During drain implantation period, tissue in growth may occur around the drain and into the holes. This could cause the drain to break upon removal. The patient's rate of healing should monitored carefully by the surgeon. If tissue in grows occur, drain's removal may need to be accompanied by surgical intervention.

All serious or life-threatening adverse events or death link to the use of these devices shall be immediately reported to the manufacturer. Any device which has been used or partially used must be destroyed in proper incinerator in accordance with the local or national regulations concerning the disposal of hospital waste.

Do not use if the packaging has been opened or damaged.

The design of this device doesn't allow its reuse.

Destroy after use. If the drain is connected to a suction drainage bottle, respect following suctions:

French Gauge	FG. 10	FG. 15	FG. 19	FG. 24	Width 7 mm et 10 mm
Recommended suction (mbar)	950	700	500	400	950

**SHELF LIFE**

Validity 5 years after the date of sterilization.

**SYMBOLS USED ON LABELLING**

	= Use by date		= Caution		= Keep dry
	= Batch code		= Do not re-sterilize		= Manufacturer
	= Catalogue number		= Do not re-use		= Keep away from sunlight
	= Authorized representative in the European Community		= Do not use if package is damage		= See instructions for use
	= date of manufacture		= Green Dot package recycling		= Sterilized using ethylene oxide
	= Not made with phthalates plasticizers		= Not made with natural latex		
	= Not made with BPA				

**CE 0123** = CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC

**Distributed by/ Distribué par:**  
PETERS SURGICAL – Immeuble AURELIUM – 1 Cours de l'Île Seguin – 92100 BOULOGNE-BILLANCOURT – FRANCE  
**TEL.:** + 33 (0)1 48 10 62 62 | **FAX:** + 33 (0)1 48 91 22 99 | **E-MAIL:** peters@peters-surgical.com  
**Manufactured by:** Degania Silicone Ltd. Degania Bet, 1513000 Israel  
**CE Representative:** Biometrix S.R.O.: Vicenzy 16, Samorin 93101 - Slovakia

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Výrobca: Degania Silicone Ltd, Degania Bet, 1513000 Israel