- (FR) APPLICATEURS DE CLIPS HEMOSTATIQUES VITALITEC®
- (EN) VITALITEC® HEMOSTATIC CLIP APPLIERS
- (DE) VITALITEC® APPLIKATOREN FÜR HÄMOSTATISCHE CLIPS
- (IT) APPLICATORI DI CLIP EMOSTATICHE VITALITEC®
- (PT) APlicadores de clipes hemostáticos VITALITEC®
- (ES) APlicadores de clips hemoestáticos VITALITEC®
- (NL) APPLICATOREN VOOR HEMOSTATISCHE CLIPS VAN VITALITEC®
- (DA) APPLIKATORER TIL VITALITEC® BLODSTANDSENDE KLEMMER
- (SV) VITALITEC® APPLIKATORER FÖR HEMOSTATISKA LIGERINGSKLÄMMOR
- (TR) VITALITEC® HEMOSTATİK KLİP APLİKATÖRLERİ
- (SK) APLIKÁTORY HEMOSTATICKÝCH SVORIEK VITALITEC®
- (RU) ГЕМОСТАТИЧЕСКИЕ КЛИП-АППЛИКАТОРЫ VITALITEC®
- (JA) VITALITEC® 止血用クリップアプライヤー
- (AR) قصبان طبية لوضع الملاقط المساعدة على تخثر الدم من فيتاليتيك®
- (ZH) VITALITEC® 止血钳夹施夹器

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# VITALITEC® HEMOSTATIC CLIP APPLIERS

## INSTRUCTIONS FOR USE

### PLEASE READ ALL INFORMATION CAREFULLY

Failure to properly follow the instructions may lead to serious surgical consequences.

Incorrect handling and care as well as misuse can lead to premature wear of surgical instruments and to serious damages. All persons using this device should be knowledgeable in the use and handling of surgical instruments and associated accessories and equipments.

### INDICATION

VITALITEC® Appliers are intended for the loading, holding and optimal closure of VITALITEC® hemostatic clips.

### CONTRAINDICATION

These appliers should not be used for anything other than their intended use.

### INITIAL USE

All instruments are carefully inspected before shipment. Due to damage that may occur during transport, the instruments should be thoroughly inspected upon receipt.

This product is delivered non sterile, it must be properly cleaned and sterilized before use.

### INSPECTION

Each instrument must be carefully inspected to assure proper functioning prior to each use with particular attention paid to the jaws and all moving parts, pins, springs, grooves, box locks.

#### Particular attentions:

- Make sure that the two jaws first close at the distal end and are not bent.
- Make sure the jaw alignment is correct to prevent scissoring effect and to avoid an improper clip closure
- Make sure the grooves are not warped and/or do not show any damage.

Inspect carefully the applier prior to each use, to identify any signs of unusual wear, cracking or corrosion. Microscopic cracks may be detected by the appearance of rust on the instrument. Should any anomaly be noted, DO NOT USE THE INSTRUMENT.

### HANDLING

Instruments should be handled and operated by staff fully familiar with their use. All surgical instruments should be handled with the greatest care, cleaned, sterilized and properly stored. This is especially important for the jaws, which are a very sensitive area.

In any case, for all decontamination, cleaning and sterilization procedures, staff should follow the instructions described below.

### INSTRUCTIONS FOR CLEANING, DISINFECTION AND STERILIZATION

All instruments must be decontaminated and thoroughly cleaned and sterilized after each use. Clean, inspect and test each instrument carefully. Sterilize all instruments before use. A proper cleaning and regular maintenance tests will extend the life of the appliers.

<b>GENERAL COMMENTS</b>	This product is delivered non sterile, it must be properly cleaned and sterilized before use. The cleaning and sterilization equipment must be validated and subjected to regular maintenance and controls. The validated parameters of the cleaning and sterilization cycles must also be controlled regularly. Cleaning and disinfecting products: make sure to follow the manufacturer's instructions in regard to the concentration, the time of contact and the temperature. Water: Preferably use deionized water with a low bacterial content (max. 100 germs/mL) and a low endotoxin rate (max. 0.25 UI of endotoxin/mL) for preparing the solutions and rinsing. Air: Use filtered air for drying.
<b>STAFF PROTECTION</b>	Staff should wear gloves when handling cleaning and disinfection products. Refer to the Technical Data Sheets for the products used.
<b>PREPARATION</b>	During cleaning, disinfection and sterilization, the applier must be opened as follows: - Disengage the middle spring. - Then spread handle rings completely for a maximum opening
<b>PRE-DECONTAMINATION</b>	The pre-decontamination procedure should be performed as soon as possible after applier's use. 1- Open the applier 2- Remove the main impurities immediately after use, with a single use tissue. 3- Place the applier in a disinfecting or enzymatic solution, with a neutral or weakly alkaline pH (pH<11), for the time and at the temperature recommended by the manufacturer (ex: Hexanios G+R at 0.5% during 15 minutes at room temperature) 4- Rinse abundantly with deionized water.
<b>MANUAL CLEANING</b>	1- Brush the whole applier with a soft brush and a neutral or weakly alkaline (pH<11) detergent until complete disappearance of visible soiling. (ex : Hexanios G+R at 0,5%, 40°C). 2- Rinse abundantly with deionized water during at least 2 minutes. 3- Check the absence of visible soiling.

<b>MANUAL DISINFECTION</b>	<p>1- Place the applier in a disinfecting solution, neutral or weakly alkaline (pH&lt;11), for the soaking time, at the concentration and temperature recommended by the manufacturer (ex : Hexanios G+R at 0,5%, 40°C during 15 minutes). Make sure the applier is completely immersed.</p> <p>2- Remove the applier and rinse 3 times during at least 5 minutes in a deionized water tank. The deionized water should be changed between each drenching.</p> <p>3- Dry for 30 minutes with hot air at 100°C.</p>
<b>AUTOMATIC CLEANING - DISINFECTION</b>	<p>Use a washer-disinfector validated according to the ISO 15883 Standard.</p> <p>Use a neutral or weakly alkaline (pH&lt;11) detergent suitable for the washer-disinfector following the recommendations of the manufacturer in regard to concentrations, time and temperature. The rinsing water must be deionized.</p> <p>1- Place the applier on a suitable support, the instruments must not be in contact with each other.</p> <p>2- Launch the cycle.</p> <p>Example of cycle :</p> <p>Washer-disinfector Miele G 7882, program Vario TD :</p> <ul style="list-style-type: none"> <li>- Pre-wash: 3 minutes with cold water</li> <li>- Wash: 5 minutes at 55°C with 30g of Neodisher MA detergent</li> <li>- Neutralization: 6 minutes at 25°C with 10mL of Neodisher Z.</li> <li>- Intermediate rinsing: 6 minutes with deionized water at 20°C</li> <li>- Thermal disinfection: 5 minutes at 90°C (A0=3000)</li> <li>- Drying: 35 minutes with hot air at 100°C.</li> </ul> <p>3- Remove the applier from the washer-disinfector.</p>
<b>CONTROL</b>	After the cleaning-disinfection operations, check that the applier has no marks of corrosion, impurities or damaged parts. A particular attention should be paid to the jaws of the applier.
<b>PACKAGE</b>	Place the product in a suitable packing for steam sterilization according to ISO 11607 standard.
<b>STERILIZATION</b>	<p>Use a fractionated vacuum sterilizer in conformity with ISO 17665 Standard.</p> <p>Sterilization cycle: 132°C for at least 2 minutes (SAL 10-6).</p> <p>In France, validated cycle: 134°C during 18 minutes.</p> <p>Maximum temperature: 138°C.</p>

**Note :** After each cleaning and before sterilization, it is recommended that the appliers, and especially the box lock and all moving parts, be treated with an approved surgical instrument lubricant.  
Make sure the lubricant is compatible with the method of sterilization chosen and please make sure to follow the manufacturer's instructions for lubricant dilution and handling.

## STORAGE

Store the products in a clean dry area with good ventilation away from areas of extreme temperature and humidity.

## SERVICE / MAINTENANCE

Never attempt to make repairs yourself. Instruments returned to Peters SURGICAL for service must have a statement certifying that each instrument has been thoroughly cleaned and sterilized.  
Service should be referred to Peters SURGICAL trained qualified staff only.

The instrument life is not specified. Designed and manufactured according to the most stringent standards, the instruments will perform for a reasonable number of years if these instructions are properly followed.

The manufacturer declines responsibility for problems resulting from improper use.

## FOR UK ONLY

**AUTOMATIC CLEANING - DISINFECTION:** use a washer-disinfector validated to HTM2030.  
**STERILIZATION:** use a fractionated vacuum sterilizer in conformity with HTM2010.

<b>MD</b>	Medical device	<b>REF</b>	Catalogue number
<b>LOT</b>	Batch code		Manufacturer
	Date of manufacture (YYYY-MM-DD)		Non-sterile
	Keep dry		Unique Device Identifier
	Caution		(For USA): CAUTION: US Federal Law restricts this device to prescription only.
	Consult instructions for use		CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC

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