Peters SURGICAL

Immeuble AURELIUM
1 Cours de l'Île Seguin
92100 BOULOGNE-BILLANCOURT
FRANCE
Tél. +33 (0)1 48 10 62 62
Fax +33 (0)1 48 91 22 99
E-mail : peters@peters-surgical.com
www.peters-surgical.com

MD

REF

AIG2120 - AIG2150

Fr

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CAUTION

Please read these instructions carefully before using PETERS SURGICAL disposable pneumoperitoneum needles. Failure to comply with the following recommendations could result in serious consequences. Laparoscopic operations must be carried out solely by doctors who have been specifically trained and have sufficient experience to use this technique. Indeed, every surgeon who undertakes an laparoscopic operation must be aware of the risks and possible complications involved with this surgical technique (notably risks to the patient and to devices involving the use of any external power source, such as electrosurgery, ultrasound, laser, gas and liquid pressure). Ensure that the patient has been informed about the way the operation will take place, in accordance with current legislation.

DESCRIPTION

PETERS SURGICAL disposable pneumoperitoneum needles have an external diameter of 2 mm. The length of PETERS SURGICAL disposable pneumoperitoneum needles is 120 mm or 150 mm, as indicated on the product label. PETERS SURGICAL disposable pneumoperitoneum needles consist of a traumatic beveled needle and a stylet of the same length with a blunt tip. The chuck is linked to a red position indicator. A Luer Lock provides connection to the gas source. A one-way valve opens and closes the gas inlet.

INTENDED USE OF PNEUMOPERITONEUM NEEDLES

The PETERS SURGICAL disposable pneumoperitoneum needle is a single-use, sterile device intended for percutaneous insertion in the

peritoneal cavity for the purpose of insufflation with carbon dioxide to establish a pneumoperitoneum prior the placement of trocars during laparoscopic procedures.
They must be used by duly qualified and experienced practitioners.

CONTRAINDICATIONS

The use of PETERS SURGICAL pneumoperitoneum needles is strictly excluded when there is a contraindication for laparoscopic surgery.

- **PETERS SURGICAL can in no case be held responsible for the assessment of any contradiction for any laparoscopic operation.**

PRECAUTIONS

- Procedure-related precautions are repeated in the directions for use section (in bold and preceded by the symbol: **!**).
- These instructions solely contain information concerning PETERS SURGICAL pneumoperitoneum needles and how to use them. They must in no way be considered as general instructions concerning techniques for the use of pneumoperitoneum needles. The surgeon alone is responsible for using PETERS SURGICAL disposable pneumoperitoneum needles in suitable fashion, bearing the recommendations in mind.
- PETERS SURGICAL pneumoperitoneum needles are designed for single use only. They must in no way and on no account be used a second time.
- The pneumoperitoneum needles are delivered in sterile packaging. They must not be resterilised, regardless of whether or not they have been used.
- The reuse or resterilisation of this device compromises its sterility and may lead to infection or cross contamination between patients. Reuse or resterilisation can also compromise the integrity and functioning of the device and may lead to an injury or pathology or to the death of the patient.
- PETERS SURGICAL shall in no case be held responsible if the rules for good surgical practice are not applied by the operator with respect to checking the compatibility with other laparoscopic instruments used during the procedures.
- When not in use, the user should pay strict attention not to injure himself when handling the needles (sharp object).
- Responsibility for the storage, opening, suitable use, manipulation and disposal of the pneumoperitoneum needles lies solely with the operator. PETERS SURGICAL shall in no case be held responsible if the rules for good surgical practice are not applied by the operator with respect to the storage, opening, suitable use, manipulation and disposal of the pneumoperitoneum needles. The manufacturer is likewise relieved of all responsibility in the case of reuse and/or resterilization.
- Do not use but destroy the sterile package together with its contents in the following situations:
 - Date marked on packaging has been exceeded
 - Packaging is severely crushed or clearly damaged.
 - Traces of humidity are visible on the packaging.
 - The insulation sheath is damaged.
- **Damaged pneumoperitoneum needles must not be used**

STORAGE CONDITIONS

- Storage temperature : avoid any prolonged exposure to low or high temperatures.
- Store in an environment that is suitable for sterile medical devices.
- During transport, apply all instructions above.

Information for the United States (USA) :

- Federal US law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- Sterilized by ethylene oxide.
- Single use only – do not reuse.

DIRECTIONS FOR USE / INSERTION PROCEDURE

1. Open the multipack box and take out the medical device.
2. Check that the packaging ensuring sterile conditions for the medical device is neither open nor damaged. If there is the slightest doubt, destroy the sterile package and its contents.
3. Check the expiry date for use marked on the package.
4. Open sterile packaging.
 - **PETERS SURGICAL pneumoperitoneum needles and their packaging must be destroyed if sterility is not guaranteed; notably in the following cases: if the date marked on the packaging has been exceeded, the packaging has been severely crushed or is clearly damaged, traces of humidity are visible on the packaging, the pack is punctured (sterile barrier).**
5. Place the PETERS SURGICAL disposable pneumoperitoneum needles on a sterile drape, taking care not to damage the needle by a blow or similar movement.
 - **Damaged pneumoperitoneum needles must not be used.**
6. Select the point of insertion. Make the incision using the usual technique used for laparoscopic surgery.
 - **Before inserting the PETERS SURGICAL disposable pneumoperitoneum needle, the operating field – the abdominal wall – must be prepared according to the usual techniques. The patient's individual anatomy must be checked before the PETERS SURGICAL disposable pneumoperitoneum needle is inserted in order to avoid any vascular lesions or injury to internal organs during insertion.**
 - **Work with the usual surgical techniques in order to protect the vessels and organs from possible damage.**
7. Insert the needle slowly, controlling the movement carefully. During insertion, especially when user holds the body of the needles, make sure that red indicator is visible at all times to monitor passage through tissue layers.
8. When the PETERS SURGICAL pneumoperitoneum needle is in place, check its position inside the abdomen.

- **The incorrect positioning of a pneumoperitoneal needle for the purposes of insufflation may lead to sub-cutaneous emphysema or embolism.**

9. Connect the insufflation tubing to the needle using the Luer Lock connection piece.
10. The Luer Lock valve has to be used in two positions :
 - a- open : to allow gas to pass through,
 - b- closed : to prevent gas to go through the pneumoperitoneum needle.
11. Use the usual safety protocol for positioning the pneumoperitoneum needle and for laparoscopic surgery techniques.
12. Continue with the procedure and remove the PETERS SURGICAL disposable pneumoperitoneum needle according to the usual laparoscopic surgery techniques.

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