DARVIN LOC 90

CE 0373

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INSTRUCTIONS FOR USE

Description:

Darvin Loc 90 suture is a sterile synthetic absorbable surgical suture device, made from poly(glicolide-co-ε-Caprolactone). The thread is available both white natural and dyed – violet (D&C Violet N°2 – color index 60725). Darvin Loc 90 sutures are available in a range of gauges sizes (USP/EP) and lengths, non-needled or attached to stainless steel needles of different size and shapes and are barbed. Can be manufactured with a button (made of same material used for the suture) to anchor the thread.

Indications:

Darvin Loc 90 sutures are intended for use in all surgical procedures, at the user's discretion.

Properties:

The progressive loss of the suture tensile strength occurs in consequence of the poly(glicolide-co-ε-Caprolactone) hydrolysis, which allow absorption and subsequent metabolization in the body.

Studies on rats have proved that the absorption of Darvin Loc 90 suture is complete within 90-120 days from implantation.

POST IMPLANTATION

1 WEEK

2 WEEKS

APPROX. TENSILE STRENGTH RESIDUAL

~50%

~20%

Precautions:

Darvin Loc 90 sutures are to be used in hospitals operating theatres and ambulatories or in any hospital environment by qualified personnel, experienced in technical and surgical procedures where the use of an absorbable suture is required. It is therefore recommended that the surgeon will consider the possible use of non absorbable sutures for the closure of tissues which may undergo distension or stretching or requiring additional support. This suture kind may be inappropriate in elderly, malnourished or weakened patients or in patients suffering from conditions which may delay wound healing process. Care should be taken when handling the suture material to avoid any bending or flattening to the thread due to the application of surgical instruments such as forceps or needleholders. Surgical sutures are disposable medical device and must be disposed of as hospital waste material.

Adverse reactions:

Darvin Loc 90 suture may be cause local inflammatory reaction.

Warnings:

- Disposable
- Do not use if package is damaged
- Use by expiry date
- Keep in cool dry place, away from any direct heat

Sterility:

- Darvin Loc 90 sutures are sterilized by Ethylene Oxide
- Sterile unless package is opened or damaged
- Do not re-sterilize

Symbols used on labelling:

= Do not reuse

= Use by Year and Month

= Manufacturer

= Sterile unless package is opened or damaged. Sterilization Method: Ethylene Oxyde

= CE mark and identification number of the Notified body. Product conforms to the essential requirements of Directives 93/42/EEC and 2007/47/EC concerning Medical Devices

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