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
Français (French) .....
English .....
Español (Spanish) .....
Deutsch (German) .....
Português (Portuguese) .....
Italiano (Italian) .....
Ελληνικά (Greek) .....
Svenska (Swedish) .....
Dansk (Danish) .....
Norsk (Norwegian) .....
Nederlands (Dutch) .....
Русский (Russian) .....
عربى (Arabic) .....
Polski (Polish) .....
Magyar (Hungarian) .....
Türkçe (Turkish) .....
中文 (Chinese) .....
Romană (Romanian) .....
Ceská (Czech) .....

APPLICATEUR POUR COLLE IFABOND® .....	2
IFABOND® GLUE APPLICATOR .....	3
APLICADOR PARA ADHESIVO IFABOND® .....	4
APPLIKATOR FÜR CHIRURGISCHEN KLEBSTOFF IFABOND® .....	5
APLICADOR PARA COLA IFABOND® .....	6
APPLICATORE PER COLLA IFABOND® .....	7
ΑΠΛΙΚΑΤΕΡ ΓΙΑ ΚΟΛΛΑ IFABOND® .....	8
MUNSTYCKE FÖR IFABOND®-LIM .....	9
APPLIKATOR TIL IFABOND® LIM .....	10
APPLIKATOR TIL IFABOND®-LIM .....	11
APPLICATOR VOOR IFABOND® LIJM .....	12
АППЛИКАТОР ДЛЯ КЛЕЯ IFABOND® .....	13
IFABOND® اپلیکاتور چسب .....	14
APLIKATOR DO KLEJU IFABOND® .....	15
ADAGOLÓ AZ IFABOND® RAGASZTÓHOZ .....	16
IFABOND® YAPIŞTIRICI İÇİN APLİKATÖR .....	17
IFABOND® 胶合剂用涂抹器 .....	18
APLICATOR PENTRU ADEZIV IFABOND® .....	19
APLIKÁTOR LEPIDLA IFABOND® .....	20

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PRODUIT STERILE A USAGE UNIQUE

**DESCRIPTION**

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STERILE SINGLE USE PRODUCT

## DESCRIPTION

IFABOND® glue applicators are class IIa medical devices. They are composed of a PVCLuer lock adapter and two hollow PVC tubes: a transparent external tube with a diameter that enables the use of the applicators with a 5-mm trocar and a green internal tube with an internal diameter of 0.7mm for standard applicators, or a transparent internal tube with an internal diameter of 0.4 mm for drop-by-drop applicators thereby providing two different delivery methods.

The applicators are available in three different lengths: 15 cm (for open surgery), 37 cm (laparoscopic surgery) and 45 cm (for bariatric surgery). The syringe provided with the glue kit is connected with the Luer lock adapter which also ensures that the cannula is securely locked.

Due to their specific design, the Drop-by-drop applicators enable perfect control of flow rate and to apply the IFABOND® glue drop by drop.

For standard applicators, IFABOND® glue can be delivered in a uniform line, thereby enforcing suture lines.

The tube is made of a malleable material and enables precise application of the IFABOND® surgical glue.

<b>MB37S</b>	STANDARD GLUE APPLICATOR - 37 CM
<b>MB15G</b>	DROP-BY-DROP GLUE APPLICATOR - 15 CM
<b>MB37G</b>	DROP-BY-DROP GLUE APPLICATOR - 37 CM
<b>MB45G</b>	DROP-BY-DROP GLUE APPLICATOR - 45 CM

## INDICATIONS

The MB devices are designed to control the application of IFABOND® surgical glue during open surgical or coeliosurgical procedures.

These devices are used to enable the adherence of prostheses or tissues and for haemostasis, sealing or suture reinforcement.

## WARNING/PRECAUTION

Compatibility studies have been conducted with IFABOND® surgical glue. It is therefore strongly advised not to use these applicators with other types of glue and to read the instructions for use of surgical glues before using them.

The applicators should only be used by experienced persons who are trained in the use of this type of device.

The use of standard applicators requires that the user control the quantity of glue delivered to avoid placing too much on the same spot.

During the procedure, the applicator tip may come in contact with blood or other organic fluids which could result in the polymerisation of the glue on the tip thereby blocking the applicator.

The applicators should not be used if the device is damaged.

Do not reuse, re-treat or re-sterilise. Reuse, re-treatment or re-sterilisation can compromise the structural integrity of the device and/or lead to a defect in the device that could cause a lesion, disease or death of the patient. Reuse, re-treatment or re-sterilisation could also result in the contamination of the device and/or cause an infection or cross infection in the patient, including notably, the transmission of infectious diseases from one patient to another. Contamination of the device can cause a lesion, disease or death of the patient.

Any used or opened device should be destroyed in compliance with local or national regulations concerning the disposal of hospital waste.

## USE

- Use sterile technique when handling the applicator.
- Mould the applicator tube to the most suitable shape.
- Fill the syringes beforehand with IFABOND® surgical glue and then lock the syringe to the applicator.
- Apply the glue directly to the desired areas by pressing gently on the piston of the syringe. The glue is then released drop by drop or in a continuous flow.
- If the applicator tip touches tissues, the polymerisation of the glue can block the lumen of the tube. The tip must be immediately cleaned with a damp compress until the polymerised glue has been removed.

## STERILISATION

MB applicators are sterilised in a peel pouch with ethylene oxide.

## PACKAGING

Sterile, in a peel pouch, in boxes of 12 single-use units  
Before opening the package please inspect it to ensure soundness.

## EXPIRY/STORAGE CONDITIONS

Do not use after the expiry date

This device should be kept in its original package in a dry place, sheltered from light and between 0°C and 40°C.

## SYMBOLS USED ON THE PACKAGE

	= Medical device
	= Do not reuse
	= Do not resterilize
	= Use by date (YYYY-MM-DD)
	= Date of manufacture (YYYY-MM-DD)
	= Manufacturer
	= Batch number
	= Do not use if package is opened or damaged
	= Keep away from sunlight
	= Keep dry
	= Upper temperature limit
	= Does not contain latex
	= Does not contain phthalates or derivatives
	= Sterilized using ethylene oxide
	= Single sterile barrier system
	= catalogue number
	= Unique Device Identifier
	= Caution
	= Consult instructions for use
	= CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC

The manufacturer should be notified of any serious adverse events, life-threatening events or events that lead to death resulting from the use of this device.

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