ACIER

CE 0459

Français	Suture chirurgicale non résorbable, monofil d'acier inoxydable
English	Non absorbable surgical suture, stainless steel
Español	Sutura quirúrgica no absorbible, acero inoxidable
Italiano	Sutura chirurgica non assorbibile in acciaio inossidabile
Deutsch	Nicht resorbierbares chirurgisches Nahtmaterial, Edelstahl
Nederlandse	Niet-absorbeerbare chirurgische hechtdraad, roestvast staal
Ελληνικο	Μη απορροφησιμο χειρουργικο ραμμα, απο ανοξειδωτο χαλυβα
Polski	Nic chirurgiczna niewchłaniałna ze stali nierdzewnej
Română	Fir chirurgical de sutură neabsorbabil, din oțel inoxidabil
Ceská	Nevstřebatelný chirurgický šicí materiál z nerezavějící oceli
Magyar	Nem felszívódó, monofil, rozsdamentes acél varróanyag
Русский	Нерассасывающийся хирургический шовный материал из нержавеющей стали
عربي	خيوط الغرز الجراحية غير القابلة للامتصاص، من الصلب غير القابل للصدأ
中文	不可吸收外科手术缝线，不锈钢

PETERS SURGICAL

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ACIER CE 0459**NON ABSORBABLE SURGICAL SUTURE, STAINLESS STEEL****DESCRIPTION**

The synthetic non-absorbable surgical sutures **ACIER** are composed of implantable stainless steel wire

ACIER sutures comply with the requirements of the United States Pharmacopoeia (USP) and the European Pharmacopoeia (EP) for non-absorbable surgical sutures, except for some diameters.

INDICATIONS

ACIER sutures are intended for use in thoracic and cardiothoracic surgery (sternum closure) and in orthopaedic surgery.

ACTIONS

ACIER sutures elicit a minimal acute inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissue.

ACIER sutures are not absorbed and there is no degradation or loss of strength.

CONTRAINDICATIONS

Even if stainless steel complies with implantable steel standard it contains nickel and it's not recommended for use on nickel sensitive subjects.

PRECAUTIONS

As with any suture, care should be taken in handling to avoid damage. Avoid thread crushing or crimping damage due to application of surgical instrument, such as needle holders or forceps.

Sutures should be handled preferably with a needle holder at a third or half of the distance between the connection to the braid and the point to avoid penetration difficulties, rupture or needle bending. Changing the bending of the needle leads to a torsion resistance loss and can cause the breakage.

The design of this device doesn't allow its reuse.

Do not resterilize. Sterile unless packaging has been opened or damaged. Destroy after use.

Users should be familiar with surgical procedures and techniques using non-absorbable sutures employing **ACIER**, as the risk of wound dehiscence may vary with the site of implantation and the suture material used.

Acceptable surgical practice must be followed with respect to drainage and closure of infected or contaminated wounds.

Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury and getting infectious illness.

CAUTION: Federal (USA) Law restricts this device to sale by or on the order of a physician.

Products should be destroyed in respect to local and national hospital wastes regulation.

ADVERSE EVENTS

To perform MRI in subjects carrying **ACIER** it is necessary to have a perfect knowledge of the technique's contraindications and artefacts associated with the presence of metallic components.

Adverse effects associated with the use of this medical device include:

- wound dehiscence,
- risk, if existing microbiological infection, that suture, as any foreign substance, accentuate infection,
- minimal acute inflammatory reaction typical of a moderate inflammatory response to a foreign body,
- transitory local irritation.

HOW SUPPLIED

The synthetic non-absorbable surgical sutures **ACIER** are available in various lengths, diameters (USP/EP) and quantities with or without surgical needles.

SHELF LIFE / PREVENTION

Recommended storage conditions: below 40°C, away from moisture and direct heat.

Do not use after expiration date.

SYMBOLS USED ON LABELLING

	= Medical device
	= Unique Device Identifier
	= catalogue number
	= Do not reuse
	= Do not resterilize
	= Use by date (YYYY-MM-DD)
	= Date of manufacture (YYYY-MM-DD)
	= Manufacturer
	= Batch number
	= Do not use if package is opened or damaged
	= Keep away from sunlight
	= Keep dry
	= Upper temperature limit
	= (for USA): CAUTION: US Federal Law restricts this device to prescription only.
	= Sterilized using ethylene oxide
	= Single sterile barrier system
	= Caution
	= Consult instructions for use
	= CE Mark. Product conforms to the essential requirements of the European Medical Devices Directive 93/42/EEC

All serious or life-threatening adverse events or deaths associated with use of this device should be reported to the manufacturer.

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