**INTRACK® Surgical Clamp**
**CE**
**INTRACK® Inserts**
**CE**<sub>...</sub>

					
Latex Free	Prescription Only	Do not use if sterile package is opened or damaged	Manufacturer	Do not resterilize	Sterilized using irradiation

**INDICATIONS FOR USE**

For Vascular

- Suitable for veins and arteries
- Performs efficiently on either diseased or normal vessels
- Provides occlusion of atherosclerotic vessels without excessive closing forces
- Minimizes intimal damage and fragmentation of atherosclerotic material
- Can be clamped over indwelling catheters

For Gastrointestinal

- Replaces bulky, rubber shod clamps
- Cushion-design of inserts enables occlusion without crushing the bowel

**DESCRIPTION**

The replaceable **INTRACK®** inserts are adaptable to any shape jaw within the clamp **INTRACK®** series.

**INSTRUCTIONS FOR USE**

**Placement and Removal of Inserts:**

- Hold the clamp in one hand and open the clamp jaws.
- Place the angled edge of the insert directly into the channel using the thumb and forefinger of your free hand. Continue pushing the insert through the channel of the jaw until the incision is firmly in position at the tip of the jaw. See Figure 1.
- Insert removal technique applies to all **INTRACK®** inserts. The inserts are disposable and must be removed before cleaning or sterilization of the clamp.
- The inserts can be removed by hand using the thumb on the gripping surface or by grasping the proximal end of the insert between the thumb and first finger and removing it from the channel of the jaw.
- Dispose of the inserts with other contaminated supplies.

**Inspection of INTRACK® Clamp:**

The clamp should be inspected carefully prior to each use and during the cleaning procedure to identify any signs of unusual wear, cracking or corrosion. Microscopic cracks may be detected by the appearance of rust on the instrument. Should any irregularity be noted, the clamp should not be used.

**Cleaning of INTRACK® Clamp:** Immediately after each procedure and before sterilization, the **INTRACK®** clamp should be thoroughly cleaned with detergent and water, and scrubbed with a soft brush to remove any remaining visible debris.

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**Sterilization INTRACK® Clamp:** The **INTRACK®** surgical clamps may be steam sterilized in instrument sets or individually in a separate peel pouch.

The following cycles, as described in ANSI/AAMI ST46: Steam Sterilization and Sterility Assurance are consistent with the sterilization of stainless steel surgical instruments.

<b>Pre-vacuum Cycle:</b> (2.5 psia pre-exposure evacuation) Exposure Temperature: 270° F Exposure Time: 4 minutes	<b>Gravity Sterilization Cycles:</b> Exposure Temperature: 250° F Exposure Time: 30 minutes Exposure Temperature: 275° F Exposure Time: 15 minutes
--	--

**DESIGN FEATURES**

Clamps:

- Supplied non-sterile (reusable)
- Channelled jaws provide ease of insertion
- Clamps available in choice of 33mm, 45mm, 66mm and 86mm jaw lengths in a variety of jaw shapes
- Comfortable inserts adapt to a variety of vessel sizes and available jaw shapes
- Inserts provided in choice of 33mm, 45mm, 66mm and 86mm, corresponding to the jaw lengths of the clamps

**HOW SUPPLIED**

The **INTRACK®** inserts are supplied sterile, two (2) per package, for single use only. **INTRACK®** surgical clamps are supplied non-sterile, one per package.

**PRECAUTION**

It is possible to inadvertently pass a needle through the jaw insert, thereby attaching the insert to the suture line. Care must be taken to avoid pinching gloved finger between the insert and the clamp channel.

**WARNINGS**

Do not resterilize and/or reuse the **INTRACK®** inserts. Resterilization and/or reuse may compromise the integrity of the device. This may lead to potential risks of failure of the device to perform as intended and/or cross-contamination associated with using inadequately cleaned and sterilized devices.

Proper disposal of contaminated devices is required, including disposal of reusable clamps that are no longer acceptable for use.

**STORAGE**

Store in an area that is dry and clean, away from heat and chemical fumes.

**CLINICAL BENEFITS**

Summary of the clinical evaluation report stating clinical benefits can be provided upon request.

**NOTICE**

If any serious injury has occurred to the patient or user of this device, notify Vitalizec or appropriate authority with the details.

**PRODUCT INFORMATION DISCLOSURE**

Vitalizec excludes all warranties, whether expressed or implied by operation of law or otherwise, including, but not limited to: Any implied warranties of MERCHANTABILITY or FITNESS, since storage and handling of this device by the user, as well as factors relating to the patient, the disease, treatment, surgical therapy, and other matters beyond Vitalizec's control directly affect this device and the results obtained from its use. Vitalizec shall not be liable for INCIDENTAL or CONSEQUENTIAL LOSS, DAMAGE, or EXPENSE directly or indirectly arising from the use of this device. Vitalizec neither assumes, nor authorizes any other person to assume for it, any other or additional LIABILITY or RESPONSIBILITY in connection with this device.

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