ENGLISH

NOVACLIP®
Surgical spring clip
INSTRUCTIONS FOR USE

LATEX	Rx Only						STERILE R
Latex Free	Prescription Only	Do not use if sterile package is opened or damaged	Manufacturer	Do not resterilize	Sterilized using irradiation		

INDICATIONS FOR USE

The NOVACLIP® Surgical Clip is indicated for use in peripheral vascular, cardiovascular and general surgery. NOVACLIP® is used for temporary occlusion of blood vessels. It is suitable for veins and arteries, provides occlusion of atherosclerotic vessels or normal vessels and may be used over indwelling catheters. NOVACLIP® may be used for temporary occlusion of the autogenous saphenous vein during coronary bypass surgery. NOVACLIP® may also be used as a suture tag.

DESCRIPTION

The NOVACLIP® Surgical Spring Clip provides secure holding and occlusion of vascular structures.

NOVACLIP® is designed for easy placement and removal when used in conjunction with a standard surgical grasper or forcep. The finger grips facilitate accurate placement and the attachment pads enable the spring clip to be placed properly per surgical technique.


INSTRUCTIONS FOR USE

The NOVACLIP® Surgical Clip is applied by pressing the finger grips between the thumb and forefinger to open the jaws, placing it over the vessel to be occluded, then slowly releasing finger pressure. Spring tension closes the clip and provides occlusive pressure. Elements of the clip are radioopaque.

An alternate method for applying the clip is to use a grasper or forcep. Compression of the grasper or forcep on the finger grips opens the jaws of the clip. When the instrument is slowly released, the spring activation will close the clip jaws.

DESIGN FEATURES

- Lightweight
- Low Profile
- Easy to apply manually or with standard surgeon's tools
- May be used for anchoring suture
- Radioopaque Components

HOW SUPPLIED

The NOVACLIP® Surgical Spring Clip is supplied sterile, one (1) per package. NOVACLIP® Surgical Spring Clip is for single use only.

WARNINGS

Do not resterilize and/or reuse the device. Resterilization and/or reuse may compromise the integrity of the device. This may lead to potential risks of failure of the device to perform as intended and/or cross-contamination associated with using inadequately cleaned and sterilized devices.

Proper disposal of contaminated devices is required.

STORAGE

Store in an area that is dry and clean, away from heat and chemical fumes.

CLINICAL BENEFITS

Summary of the clinical evaluation report stating clinical benefits can be provided upon request.

NOTICE

If any serious injury has occurred to the patient or user of this device, notify Vitalitec or appropriate authority with the details.

PRODUCT INFORMATION DISCLOSURE

Vitalitec excludes all warranties, whether expressed or implied by operation of law or otherwise, including, but not limited to: Any implied warranties of MERCHANTABILITY or FITNESS, since storage and handling of this device by the user, as well as factors relating to the patient, the diagnosis, treatment, surgical therapy, and other matters beyond Vitalitec's control directly affect this device and the results obtained from its use. Vitalitec shall not be liable for INCIDENTAL or CONSEQUENTIAL LOSS, DAMAGE, or EXPENSE directly or indirectly arising from the use of this device. Vitalitec neither assumes, nor authorizes any other person to assume for it, any other or additional LIABILITY or RESPONSIBILITY in connection with this device.

■ VITALITEC INTERNATIONAL, INC.

d/b/a Peters Surgical, USA

Plymouth, MA 02360 USA

Phone: 508-747-6033

Customer Service: 800-805-6590

www.peters-surgical.com

Logo and NOVACLIP® are registered trademarks of Vitalitec International, Inc.

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