INSTRUKTIONER FOR USE: Suture PassOR™ Suture Passing Puncture Closure Device

DESCRIPTION
The Suture PassOR™ is available in two configurations: RSG, RSGXL

INDICATIONS FOR USE / CLINICAL BENEFITS

The Suture PassOR™ is indicated to be used during endoscopic and laparoscopic surgery to facilitate the placement of sutures for secure closure of trocar sites as determined by a licensed physician familiar with the possible side effects, typical finding, limitations, indications or contraindications of performing such procedure.

The pre-assembled device minimizes complications and costs associated with post-operative port site hernia, by facilitating fast and complete peritoneal/fascial closure.

CONTRAINDICATIONS

This device is not intended for use except as indicated, and when laparoscopic and other minimally invasive surgical techniques are contraindicated.

CAUTIONS / WARNINGS

- This device should be used only by surgeons trained in proper techniques for laparoscopic surgery and other minimally invasive procedures.
- Device is sterile if the package is dry, unopened, and undamaged. DO NOT use if the pouch is damaged or the packaged seal is broken.
- Verify compatibility. Inspect the endoscope/laparoscope and trocar for sharp edges that could interfere with the components of the Suture PassOR™.
- Protect the Suture PassOR™ with protective sheath when not in use. Being a needle, it can cause injuries.
- DO NOT use in surgical procedures where the position of the Suture PassOR™cannot be clearly determined.
- DO NOT attempt to use the Suture PassOR™ if the plunger is not functioning properly.
- For single patient use only. The design of this device may not perform as intended by the manufacturer if it is reused. The manufacturer cannot guarantee the performance, safety and reliability of a reprocessed device.


INSTRUKTIONER FOR USE

The Suture PassOR™ has specific indications for use as noted above. In that context, the surgeon is best advised to use a method which his/her own practice and discretion dictates to be best for the patient, consistent with the indications and contraindications outlined above. The following instructions are recommended for the proper function of the Suture PassOR™. This is not a reference for trocar closure techniques.

Open sterile pouch, remove the Suture PassOR™ and place onto sterile field.

- Acquire a suture at least 24 inches in length. Remove the needle guard from the Suture Pass-OR™. The arrow indicator mark is aligned with and designates the beveled side of the needle. Press the plunger on the proximal end of the device to advance the grasper arms past the distal end of the needle. Place the midpoint of the suture between the grasping arms. Release the plunger to retract the suture into the needle. Pull the two lengths of suture back along the needle towards the handle end of the device and hold in place with gentle tension.
- Advance the Suture PassOR™ needle through the subcutaneous tissue and fascia. Keep the needlepoint in direct visualization after it has entered the abdominal cavity. Briefly press the plunger to release the suture. Withdraw the device.
- Position and then push the needle through the opposite side of the subcutaneous tissue and fascia. Keep the needlepoint in direct visualization after it has entered the abdominal cavity. Press the plunger to advance the grasping arms. Grasp the previously placed suture loop with the grasping arms. Release the plunger to retract this suture into the needle. Hold one of the two ends of the suture to prevent movement. Withdraw the needle to pull this suture through the tis- sue. Continue to retract this suture loop until the other end of the suture exits the grasping arms.
- Using manual traction on the suture withdraw any excess length that may have remained in the abdominal cavity. Use standard technique to tie the suture.
- Upon completion of the procedure, carefully place the used needle in a sharps bio-hazard container.

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

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