



Ifabond®  
Innovative technology for mesh  
fixation and hemostasis.

Our **surgical glue**  
solutions.



# Ifabond® surgical glue.

## Designed for surgeon.

- Immediate polymerization with visible process : glue whitening color.
- Adhesive effect after only 30 seconds<sup>1</sup>.
- Different sizes of appliers with curvable distal tip.
- No rippling effect during the mesh fixation (vs tension exerted by a thread).

## Optimized for OR staff.


- Ready to use product: no mix nor preparation of component is requested.
- All in one kit: glue vial, syringe and extraction needle included.
- 3 volumes of glue available to match with surgery needs (waste prevention).
- Significant reduction of operating time<sup>1\*</sup>.

## Secured for patient.

- Atraumatic glue, designed to prevent any transfixion during the mesh fixation.
- Resorption<sup>2</sup>: partial from 3 to 6 months, total from 6 to 12 months.
- Flexible glue<sup>3</sup> after its polymerization (N-Hexyl) absorbing potential stress.
- Slow degradation without emission of toxic products.
- Post operative pains significant reduction<sup>1\*</sup>.




## Applicators available in different dimensions according to the surgical approach.




Distal tip of each applicator can be bent for a more precise Ifabond® delivery.

## A few gluing drops are enough.

In visceral and digestive surgery, Ifabond® is used for meshes fixation whatever the approach.



In open surgery.



In laparoscopic surgery.

Gluing dots fixation.



## Indications and surgical procedures.

### MESH FIXATION

### HEMOSTASIS LYMPHOSTASIS



**DIGESTIVE AND VISCERAL SURGERY**

**INGUINAL HERNIA REPAIR**

Laparoscopic surgery, Laparotomy



**URO-GYNECOLOGICAL SURGERY**

**SACROCOLPOPEXY**

Laparoscopic surgery, Laparotomy



**BARIATRIC SURGERY**




**SLEEVE - BYPASS**



**SENOLOGIC SURGERY**

**MASTECTOMY**

**AXILLARY CURAGE**



## CLINICAL STUDIES.

### ADHESION

Inguinal Hernia: Mesh fixation.

<sup>1</sup>Moreno-Egea et al. «Is It Possible to Eliminate Sutures in Open (Lichtenstein Technique) and Laparoscopic (Totally Extraperitoneal Endoscopic) Inguinal Hernia Repair? A Randomized Controlled Trial With Tissue Adhesive (n-Hexyl-Cyanoacrylate)» Surg Innov. January 7, 2014 - Randomized single-blind monocentric trial. 208 patients randomized between two groups (Open & Laparoscopic approach) treated with Prolene sutures (control) vs Ifabond (experimental). To evaluate the effectiveness of Ifabond to treat inguinal hernias in term of pain, recurrence, operating time, analgesic consumption & morbidities. Study shows that the use of Ifabond : significantly reduced pain, and analgesics consumption, both via the open and laparoscopic approach ( $p < 0.001$ ) ; did not change the recurrence rate in either of the approaches after 1 year ; no morbidity associated with the use of Ifabond existed ; significantly reduced the mean of surgical time (12 minutes in open surgery, 13 minutes in laparoscopic surgery ;  $p < 0.001$ ). The author concludes that substituting sutures with Ifabond in open or laparoscopic inguinal hernioplasty is safe with less postoperative pain and the same possibilities of recurrence.

Sacrocolpopexy: Meshes fixation.

Dr G. Lambin, 2017, J Minim Invasive Gynecol 24(1): 41-47 «Attachment in Laparoscopic Sacrocolpopexy: A Prospective Multicenter Pilot Study».  
Dr JP. Estrade et al, 2015, «Laparoscopic sacrocolpopexy with a vaginal prosthetic adhesive» Gynecol Obstet Fertil 43(6): 419-423.

Dr RK Silveira, 2017, «Comparative study of safety and efficacy of synthetic surgical glue for mesh fixation in ventral rectopexy», Surg Endosc.

### HEMOSTASIS

**BARIATRIC SURGERY (SLEEVE GASTRECTOMY).**

Dr G. Mercier MD, PhD et al, 2017, Journal of Evaluation in Clinical Practice. «Surgical glue in laparoscopic sleeve gastrectomy: An initial experience and cost effectiveness analysis».

### LYMPHOSTASIS

Axillary curage patient with breast cancer.

Conte, M. R. Payan et al. (2013) «Synthetic glue application in breast cancer surgery to prevent lymphocele: any interest ?» Randomised multicentric study. 9th day Daniel DARGENT gynecology surgery. Lyon.

### Preclinical Data.

<sup>2</sup> Rapport NAMSA-2013/14 : « Surgical glue local tissue effects and degradation evaluation ».

<sup>3</sup> MEMO RD-11-001/2011-Std NF EN ISO 10993 : « IFABOND® properties (purity, elasticity, polymerization temperature) ».

<sup>4</sup> BIBLIO-IB-02 ADHESION -Février 2014 : « IFABOND® glue adhesion in vitro evaluation ».



## References.

Code	Volumes/Lenght	Description	Qty/box
<b>IB05</b>	<b>0,5 ml</b>	N-hexyl cyanoacrylate glue + 2.5ml Luer-lock syringe + 18G puncture needle	6
<b>IB</b>	<b>1 ml</b>		
<b>IB+</b>	<b>1,5 ml</b>		
<b>MB15G</b>	<b>15 cm</b>	Drop-by-drop applicator for Ifabond® glue	12
<b>MB37G</b>	<b>37 cm</b>		
<b>MB45G</b>	<b>45 cm</b>		

## Indications.

Ifabond® glue - Class III Medical Device - CE0459 - Manufacturer Peters Surgical  
Ifabond® surgical glue is intended to be used during surgical procedures, in open surgery and laparoscopic surgery, for its adhesive and hemostatic sealant actions. These surgeries for the adult population include digestive and visceral surgery for the treatment of hernia, bariatric surgery for sleeve-bypass, urological and gynecological surgery with sacrocolpopexy and breast surgery for mastectomy and axillary curage. The glue can be used in the pediatric population for circumcision.

### Read the instructions carefully before using the products.

Presentation to Peters Surgical employees and distributors as well as healthcare professionals.

### Peters Surgical Headquarters & French Affiliate

Immeuble AURELIUM, 1 cours de l'Île Seguin,  
92100 Boulogne-Billancourt, France  
+33 (0)1 48 10 62 62  
peters@peters-surgical.com



[peters-surgical.com](http://peters-surgical.com)